

\$5 Cash Back Offer

Thank you for purchasing the Milton 2 in 1 Microwave and Cold Water Steriliser. To reward you for your purchase we are offering a \$5 cash back offer which may be redeemed by one of two methods.

Your offer must be activated within 90 days of purchase to be eligible for a refund. For full Terms & Conditions please refer to:

www.miltonbaby.com.au/cashbackoffer

Please allow 30 days from receipt for processing of your offer.

ON-LINE REDEMPTION:

1. Apply on line by visiting www.miltonbaby.com.au/cashbackoffer and complete the necessary information.
2. Don't forget to include a scanned copy or photograph of your original receipt as proof of purchase (ensure date of purchase is visible on receipt).
3. Once verified your cash-back is delivered in the form of an EFT and successful claimants should allow 30 days from when their claim has been processed and approved to receive the money in their bank account.

REDEMPTION BY MAIL:

1. Complete the form below with all relevant information (please use capital letters for legibility).
2. Include either your original receipt or a copy of your original receipt as proof of purchase.
3. Once verified your cash-back is delivered in the form of an Electronic Funds transfer (EFT) into your nominated bank account. Successful claimants should allow 30 days from when their claim has been processed and approved to receipt of payment.
4. Mail completed form and receipt to:
MILTON 2 IN 1 STARTER KIT CASH BACK OFFER
PO BOX 630, BRAESIDE, VICTORIA, 3195.

MILTON \$5 CASH BACK OFFER

Title: _____ First Name: _____

Surname: _____

Street: _____

Suburb: _____

State: _____ Post Code: _____

Daytime Telephone: STD () _____

Email Address: _____

Your \$5 cash back will be deposited into your nominated Bank Account (Australian only) via Electronic Funds Transfer (EFT). Please complete the details below so that payment can be processed:

Account Name: _____

BSB: _____ -

Account Number: _____

- Proof of purchase included
 Date of purchase visible on receipt
 Please tick if you would like to hear more about Milton, including new products and exclusive offers. To see how Milton collects and uses your personal information, including marketing via emails and online targeted advertising, read our privacy policy at:
<http://miltonbaby.com.au/privacy-policy>

Also available from Milton:



Milton Solo Microwave or cold water Travel Steriliser

A handy, compact, unique safety vent technology single bottle steriliser that's perfect for travel, overnight stays or when you're away from home. Use with cold water or in the microwave - it's an essential for your baby bag!



Milton Anti-bacterial Tablets

Works in just 15 minutes, killing germs and bacteria. Ideal for home use or travel. No need to rinse. Can be used for home cleaning.



Milton Anti-bacterial Solution

Double concentration. Use less with the same germ kill. Hospital grade disinfectant. Multi-purpose formula, works in just 15 minutes. Can be used for baby accessories and general cleaning. The nation's favourite, used and trusted since 1916. No rinsing required.



Milton Baby Bottle Cleaner

Clinically proven to remove milk residue. Fragrance-free, colouring-free, biodegradable formula. Developed to clean all baby feeding equipment and accessories. 100% plant-based cleaning agents.



Milton Anti-bacterial Surface Spray

Cleans and disinfects all surfaces in the home, is ideal for nursery and kitchen surfaces. Leaves a lovely fresh fragrance.



Milton Anti-bacterial surface wipes

Ideal at home or on the move for quick clean-ups. 100% plant-based active ingredient. No need to rinse. Suitable for all baby equipment and accessories.

Milton Australia PTY LTD, 1/575 Darling Street, Rozelle, NSW 2039.

Call Australia 1800 506 750 or New Zealand 0800 555 895.

miltonbaby.com.au

All rights reserved. Trademarks owned by Milton International. Product colours/design may vary.
11/2018



Milton®
TRUSTED SINCE 1916

2 in 1 Combi Steriliser Starter Kit

microwave or cold water

**THANK
YOU FOR
CHOOSING
MILTON!**

FITS MOST SIZE BOTTLES
AND IN MOST MICROWAVES


IMPORTANT INSTRUCTIONS, PLEASE READ

Congratulations on choosing your new Milton Combi Steriliser with unique safety vent technology for better, faster steam sterilisation. Before using it for the first time, please take a little time to read this leaflet really carefully. Then keep it somewhere to hand, as it contains important information that you might need to refer to again.

If you do mislay the leaflet, don't worry, you can download a copy at:
miltonbaby.com.au

Before you first use your new Milton Combi, familiarise yourself with its layout and parts. Then wash all the Combi parts in warm soapy water and rinse in cold water.

PARTS GUIDE
A Pop-up safety vent
B Screw lid with submerger and seal
C Handle
D Base
E Fill line (5L) for cold water sterilisation



2 ways to sterilise choose the most convenient method for you!

1 The Cold Water Method

1 CLEAN

- Clean bottles, teats and accessories thoroughly with warm soapy water.
- Remove any milk residue. Rinse with cold water.

2 PREPARE THE MILTON SOLUTION

- Fill your Milton Combi with cold water to the 5L mark.
- Add 2 standard Milton Anti-bacterial Tablets or 30 ml (1 capful) of Milton Anti-bacterial Solution.

3 ADD ITEMS TO THE SOLUTION


- Wash your hands.
- Fully submerge items in the solution to remove any trapped air bubbles.
- When you screw on the Combi lid the integral submerger will then ensure that all items are kept in the solution.
- Screw on lid. Press the vent down if it is raised.
- Items are ready to use after **15 minutes**.
- Remove the items you need and shake off excess water – **THERE IS NO NEED TO RINSE**.
- Replace the solution every 24 hours.

IMPORTANT: Your Milton Combi is specially designed for use solely with Milton Anti-bacterial Tablets or Solution. Other sizes or makes of tablet/solution could damage your Milton Combi and affect its performance. On the product packaging of the Tablets, you may notice that you need to use 2 tablets for 4L of water. For your Combi, you only need to use 2 tablets for 5L of water. The disinfecting solution will still be effective.

WARNINGS: Avoid contact with base metal, silverware, silks, woollen or coloured fabrics (may bleach).

How to get the best from your Combi Steriliser:

- Check the power of your microwave before you begin using it.
- Dry your Milton Combi well after each use to prevent a build-up of limescale.
- For optimum effectiveness descale once a month.
- The individual parts of the Milton Combi can be cleaned in the dishwasher.
- To remove the submerger from the lid for cleaning, insert fingers into the holes as shown and twist anti-clockwise. Reverse the procedure to re-attach to the lid.
- It is only necessary to use the submerger when using the cold water sterilising method, although for convenience it can still be left attached to the lid when sterilising in a microwave.



2 The Microwave Method

1 CLEAN

- Clean bottles, teats and accessories thoroughly with warm soapy water.
- Remove any milk residue. Rinse with cold water.

2 FILL

- Pour 200ml of water in your Milton Combi base. Then load items upside down.

IMPORTANT: Your Milton Combi will not sterilise properly unless the correct amount of water is added. You can use a baby bottle to measure exactly 200ml of water.


3 LOAD AND CLOSE

- Load bottles with the opening face down.
- Screw on the lid. Press the safety vent down if raised. When steam pressure builds up inside Combi the safety vent will pop out automatically.

4 MICROWAVE

- 2 mins at 1100-1850W
4 mins at 850-1000W
6 mins at 500-800W

2 mins



5 COOL AND USE

- For safer handling, leave the Milton Combi to cool inside the microwave for at least 3 minutes after the microwave cycle has finished. Remove the Milton Combi using oven gloves to prevent scalding.
- Wash your hands well before handling sterilised items. Unscrew the lid and take the contents out using tongs.
- If you don't want to use your bottles straight away after microwaving, press the pop-up safety vent back down to keep the bottles inside the Milton Combi sterile for longer.

Did you know?

If the Milton Combi is left unopened after microwaving and the safety vent is closed, items will remain sterile in the Milton Combi for up to 24 hours.

Did you know?

With the Cold Water method of sterilisation, items can be left in the solution until you need them and will stay sterile in the solution for 24 hours. Items can be added and removed from the same solution throughout the day.

WARNINGS FOR USE IN THE MICROWAVE:

- Check the power of your microwave before you begin using it.
- Check that the accessories you want to sterilise can be microwaved.
- Always keep the microwave turntable clean and dry before use. For smaller microwaves, the turntable may need to be removed to allow the Milton Combi to fit in.
- Always ensure that your Milton Combi is clean with no old water residue before use in the microwave.
- The Milton Combi contains boiling water when in use. Handle with care and keep out of the reach of children.
- Always wait 3 minutes before removing the Milton Combi from the microwave to let it cool down. Ensure that you take care, as it will contain hot water and steam. Always place it on a hard, heat-resistant surface. Use oven gloves and once the lid is unscrewed, open it away from you to allow the steam to escape.
- Do not use Milton Anti-bacterial Tablets or Solution inside the Milton Combi when placed inside your microwave.
- Do not use metallic items in your Milton Combi as this could damage your microwave.

THIS PRODUCT MAY BREAK IF DROPPED WHILE FULL OF WATER

Ask the experts:

Why should I sterilise?

During your little one's first 12 months, their immune system is still developing so they need extra protection from germs that may cause gastroenteritis and other nasty bugs. To do this, you need to make sure that all feeding equipment is clean and properly sterilised.

How long should I sterilise for?

Right up until your baby's first birthday

What should I sterilise?

Basically everything that can go into your baby's mouth! This includes bottles, teats, soothers, plastic feeding spoons, weaning bowls, teething rings, toys as well as your breastfeeding accessories such as breastpump parts and milk storage containers.

Did you know?

At 6 months old, your baby's immune system is at its weakest and even at 1 year they will only have around 17% of adult immunity levels against the germs that cause diarrhoea and vomiting.

For more information, take a look at our website at:
miltonbaby.com.au